



PERFECTCOLL™

OPERATING MANUAL

PC SKIN TESTER

Measurement of skin moisture





Contents



| | |
|------------------------------------|----|
| System for measuring skin moisture | 4 |
| Design and technical data | 5 |
| Recommended cosmetics | 6 |
| Application | 8 |
| Care and maintenance | 10 |

System for measuring skin moisture

The skin should maintain the proper level of moisture in order to work properly. In order to maintain the skin in the best health, it is good to have the ability to easily and frequently check the level of its moisturising.

Let us introduce the innovative product, **PC SKIN TESTER**.

PC SKIN TESTER is an easy way to check the moisture content of our skin.

We wish you satisfaction
with our products

The Management Board, Scientific Committee and Employees of Perfect Coll Sp. z o.o.



Design and Technical data PC SKIN TESTER

TECHNICAL DATA

DEVICE:

| | |
|---------------------|------------------|
| Name of the product | PC SKIN TESTER |
| Dimensions | 131 x 27 x 24 mm |
| Weight | 69 g |
| Moisture | < 70 % |

POWER SUPPLY:

| | |
|---------|---------|
| Battery | 2 x AAA |
|---------|---------|

MADE IN PRC FOR PERFECT COLL



Tester cover

Display of measurement

Reset

On/ Off



Recommended cosmetics

Use **PC SKIN TESTER** only with the original cosmetics of PERFECT COLL



1

Before treatment

PC DETOX EXTRACT FROM NORWEGIAN WILD CORNFLOWER

The cleanest and most modern, all-natural extract of the flowers of Centaurea harvested from organic farming in Scandinavia, enriched with natural enzymes from Organic Kelp (organic algae). It gently cleans and soothes the sensitive skin of face, neck and the whole body. It does not contain parabens, chemical preservatives, and drying skin alcohol.

SYMBIOSIS OF SCIENCE WITH NATURE

PERFECT COLL cosmetics have unique, innovative formula developed in biotechnology, environmental laboratories in the Polish medical universities.

We are keen on nature and ecology, so the cosmetics with the PC sign are healthy for your skin.

PERFECT COLL cosmetics do not contain parabens, chemical preservatives, petroleum (paraffin), mineral oils, alcohols, and glycols drying out the skin, as well as irritating foaming agent Sodium Laureth Sulfate and other chemical substances harmful to the skin.

Refer to "Perfect Coll Products Catalogue"



2

After treatment

PC PREBIOTIC SOOTHING CREME

Natural and eco-friendly, innovative soothing cream, which due to unique composition has remarkable properties to support the work of damaged skin cells. Its components are responsible for regeneration of even the most delicate skin, and enzymes from papaya fruit provide the effect of light glowing. It is the only product on the market that contains prebiotic enzymes, which provide very high absorbency.



Application

PC SKIN

TESTER



SIMPLE STEPS TO USE PC SKIN TESTER:

- Remove the cover from the tester.
- Press the On/Off switch.
- When you hear the signal, hold the head to the skin for approx. 10 seconds.
- The device will beep signalling the completion of taking the measurement.
- The measured moisture content will be displayed on LCD.

Further instructions

- If LCD shows "LO", turn off the device, remove the batteries and insert them again.



- If the result of the measurement showed "UUU" on LCD, it will mean an error. Please measure again by pressing the Reset button.



Do not measure the moisture content on wet face.



Frequency of treatment:
No restrictions



Treatment time:
10 seconds



Care and maintenance

- Keep the device clean and dry.
- Do not heat the device.
- Replace the batteries after a month of use.
- Clean the head immediately after and before use (spray cosmetic tissue or pad with 2% salicylic alcohol).



Keep the device clean and dry.







Centrum Operacyjno-Rozwojowe Perfect Coll Sp. z o.o.

Pomorski Park Naukowo-Technologiczny (lok. A 210)
Al. Zwycięstwa 96/98, 81-451 Gdynia



kontakt@perfect-coll.pl



+48 58 767 42 24



www.perfect-coll.pl

